


# Project QuantiCellPhe- Overview

Ziqing Ye

Trcek Lab @ Johns Hopkins University


# Overview

- Goal: quantify the number of pole cells that has sufficient granule expression (?)
  - Mechanism:
    - Auto Gradient Normalization ©
    - Smart Thresholding ©
    - Connected component analysis & record nuclei coordinates
    - 3D Cell Count (and phenotype identification) ©
- 
- 2D Nuclei  
identification


# Auto Gradient Normalize

Sample image: Embryo\_3\_w2iSIM-405-DAPI\_cmle-0034

---



Original




Gradient Normalized


# Auto Gradient Normalize (Threshold result)

Sample image: Embryo\_3\_w2iSIM-405-DAPI\_cmle-0034, threshold = 55

---




Original



Gradient Normalized

# Auto Gradient Normalize – Mechanism


- Step 1: calculate the average pixel intensity of the center strip
- Step 2: adjust the intensity of each horizontal strip accordingly
- Same for vertical adjustments



*I wish it was this simple  
but nope*


# Determining the optimal threshold

- The optimal threshold for selected images across the spectrum of z values (slice 1, 34, 58, 84, 125, 128, 138) were determined through hand-tuning
- These, along with the “center average intensity” data calculated by the Auto Gradient Normalize algorithm, were plotted in Excel and underwent linear fit.
  - We got the optimal threshold as a function of center average intensity




# Smart Threshold – introduction

- Having a single threshold value for an entire slice did not work even with Gradient Normalize. (Surprise!!!!!!)




Original image



after threshold (single value, 55)

# Smart Threshold – Mechanism


- Start from a lower threshold
- Identify continuous regions that are not properly segmented under the current threshold (based on the region size), then lower the pixel value of those regions, then re-iterate until all parts are properly segmented





# Smart Threshold – result

- It is not perfect, but it performs significantly better at the tip region where the phenotypic cells are located (compared to single value thresholding).





Original single-value threshold




Smart Threshold (max iteration = 30)

# 3D Nucleus Detection

- The stack is segmented one z-slice at a time
- Each distinct bright spot on the segmented image is identified as a temporary “cell”: an object from the custom Cell class, which has x coordinate, y coordinate, z coordinate, radius, thickness, and phenotype parameters
- Multiple lists were created to store different kinds of Cell objects.
  - currentCells (all the “cells” on a single z slice, temporary)
  - allCells (all of the “permanent” cells)
  - oldCells (all of the permanent cells that were present on the previous z-slice)




# 3D Nucleus Detection



# Phenotype Check

- The phenotype detection area (white circle)'s radius is arbitrarily set to be twice the radius of the cell nucleus (black circle).
- The pixel intensity within the circles are obtained; if more than a certain percentage of the pixels have a pixel intensity of higher than a certain threshold, the cell object in that same location is marked as phenotypic




# Acknowledgements

- All image used = copyright @ Trcek Lab @ Johns Hopkins University
- Referenced a glare detection algorithm by amittn on github for spot detection syntaxes, as well as various online library documentations
- Multiple python libraries were used, including opencv, the tifffile, xlsxwriter, numpy, etc.



Tatjana Trcek (PI)



Other members of Trcek Lab for providing imaging data & testing the program